



CACI NON SURGICAL FACIAL TONING & BODY SCULPTING

CACI QUANTUM



CACI QUANTUM

Non Surgical Facial Toning & Body Sculpting System

The result of over 20 years research, the award winning CACI Facial Toning treatment has been voted* the most effective anti-ageing treatment available.

The CACI Quantum is a multifunctional system combining CACI Non Surgical Facial Toning with Body Sculpting. The system's ability to target face and body treatment areas at the same time enables you to increase your revenue without the need to increase treatment time. The CACI Quantum redefines facial muscles, reduces lines and wrinkles and tones sagging body muscles.

FARADIC

Faradic technology uses a higher intensity output than microcurrent stimulation, and so is ideal for firming and toning the larger body muscles. The current is delivered in a series of short pulses through conductive pads that are positioned over specific muscle groups. As the current passes through the pads it stimulates involuntary muscle contractions thereby passively exercising and toning the muscles.

Faradic technology enables clients to experience a total body workout without even breaking a sweat. A series of treatments will redefine body contours and provide measurable inch loss.




A Complete
Non Surgical Face and Body
Treatment Solution



FARADIC BENEFITS

- Firms and tones body muscles
- Measurable inch loss
- Stimulates blood circulation
- Passively exercises the muscles
- Tones Abs: a typical 20 minute stomach treatment is the equivalent of doing 300 perfect sit ups



MICROCURRENT

Originally developed for treating facial palsy, microcurrent therapy is now widely used as an aesthetic anti-ageing treatment. True microcurrent uses a current with an intensity of less than one millionth of an amp (<1mA). A higher intensity current will cause contractions of the facial muscles which can worsen lines and wrinkles, whereas microcurrent stimulation does not contract muscles, but instead works by a process called muscle re-education to achieve results.

Due to its dramatic lifting and toning effect, the CACI microcurrent treatment has become known as the "Non Surgical Face Lift".

HYDRATONE

Hydratone combines active microcurrent rollers with Hydro Mask, a unique, electrically conductive silicon gel mask that has been infused with powerful hydrating properties.

The rollers gently massage the face causing the mask to become electrically charged so that the whole face is bathed with rejuvenating energy.

The current delivered by the rollers recharges and replenishes the bio-electrical energy within the facial muscles, restoring tone and firmness.



MICROCURRENT BENEFITS

- Lifts and tones facial muscles
- Reduces the appearance of lines and wrinkles
- Increases blood and lymph circulation
- Enhances the penetration of skin care formulations
- Increases the production of collagen and elastin
- Increases protein synthesis and cell membrane transport
- Increases mitochondria activity and adenosine triphosphate (ATP) production


INCREASE YOUR REVENUE WITH CACI QUANTUM

When you buy the CACI Quantum, you buy more than just a treatment system, you become part of a winning team.

Share in the success of the CACI brand name and increase your revenue.

*Based on a 6 day week and 1hr treatment at £50

WEEKLY REVENUE



ANNUAL REVENUE





BODY TONING

The CACI Quantum features specific programs for slimming, toning, body contouring, breast enhancement, buttock and thigh firming.

The system's lymphatic drainage mode utilises sophisticated bio-feedback technology to deliver electrical impulses in sync with each client's heartbeat.

Unique contour discs are used in combination with faradic stimulation to specifically target difficult to tone areas like the bust, thighs and buttocks.



72 PRESET PROGRAMS

TREATMENTS

- Facial Toning ✓
- Lines / Wrinkles ✓
- Eye Bags ✓
- Skin Hydration ✓
- Neck Treatments ✓
- Body Toning ✓
- Cellulite (with ECM)* ✓
- Slimming ✓
- Hand Treatments ✓

* Optional Attachment

ELECTRO BUDS



CACI systems utilise specially developed serum filled cotton buds to enhance treatment results. The serum contains Hexapeptides and powerful anti-ageing ingredients to provide an instant face firming effect. Electro buds are also perfect for retail.

ELECTRO CELLULITE MASSAGER (ECM)



The ECM is an optional attachment that can be used in conjunction with all CACI systems that incorporate microcurrent technology. The ECM attachment combines the therapeutic benefits of massage with the effectiveness of microcurrent therapy.

Conductive rollers are electrically charged with microcurrent to gently knead tired shoulders, soothe aching limbs and firm sagging muscles. The stimulatory action of the ECM tones muscles and reduces inches by improving lymph drainage and breaking down fatty deposits.

“ The marketing and press coverage of your product is impressive. CACI is now enjoying an even higher profile and we are proud to be on board with you. ”

Anita Wilson
Beautique, Worcester



THE CACI DIFFERENCE

All CACI systems deliver microcurrent energy using ultra low frequencies and the power of a unique digitally simulated tidal wave.

This special digital waveform optimises results and comfort by gradually increasing in intensity to overcome skin resistance (Ohm's Law). It then rises to a peak in order to deliver the full power of the microcurrent energy deep into the muscle.



Touch Screen Interface

TECHNICAL SPECIFICATION

Dimensions(cm) 45W x 42D x 15H
 Weight 9Kg
 Touch Screen 10.4in

Supply Voltage 60-240 volts
 Output Voltage: 30 volts max (loaded)
 Consumption 110 vA
 Safety Class 1 type BF in accordance with IEC 601 and IEC 601-2-10
 Output Channels: 16 independent
 Output Range: 20 uA to 20 mA

MICROCURRENT

Polarity: positive, negative & alternating
 Range: 20 uA to 640 uA
 Frequency: 0.1Hz - 800Hz
 Waveshapes: 8
 Mode: Automatic Sequential

INSIDE THE BOX

ATTACHMENTS

- 1x Pair Micro Probes
- 1x Mains Lead
- 1x Slimtone Leads (Pack of 12)
- 1x T-Bar
- 1x Pair Straight Bars
- 2x Contour Heads
- 1x Lymph Cardiac Sensor
- 1x Channel Outlet Tester
- 1x Pair Hydratone Face Rollers
- 1x Pair Double Tip Probe Inserts

CONSUMABLES

- 8x Self-Adhesive Pads -Small (Set of 4)
- 16x Self-Adhesive Pads -Large (Set of 4)
- 1x Collagen Revital Serum (30g)
- 1x Essential Cleansing Milk (300ml)
- 1x Eye Wrinkle Gel (100g)
- 1x Firming Gel Concentrate (300ml)
- 1x Hydro Moisturising Cream (300ml)
- 1x Revitalising Toning Fluid (300ml)
- 1x Pack Electro-Buds (Per 100)
- 20x Hydro Masks

PROMO/MARKETING MATERIAL

- 1x A1 Poster
- 1x DL Client Leaflets (per pack)

QNTBRO